



Flow Chemistry Specialist

LTF Representative in India

AMAR EQUIPMENTS PVT. LTD.

Serving Industry since 1974.

6, Parmar Industrial Estate, Bail Bazaar, Kale Marg, Kurla (W), Mumbai - 400 070, India.

Tel: +91-22-64562888 / 25129273 / 25111540 / 25143142 / 25140169 • Fax: +91-22-25109565

E-mail: sales@amarequip.com / export@amarequip.com • Website: www.amarequip.com


Little Things
FACTORY

Plan Optik Group

Syringe pump Mr. Q

Continuous flow

from 80 μ l/h bis 3,5 liter/h (depends on syringe diameter)



Pressure sensor optional

Max. pressure 145 psi

Fluids in contact with Hastelloy and PTFE

Syringe diameter from 0,7 mm to 23 mm

Self-priming

Dispense und multidispense

pipet und multipipet

Auto restart

PC-Software for Microsoft,

Apple and Linux

Suitable for Flow Chemistry Starter Set